



**Biomatlante**  
Biologics Solutions

eXpert in bone regeneration

# EZ CURE™

## Collagen Membrane

Adaptable

Flexible

Crosslinked for long-lasting barrier





### Barrier Function lasts for 3-6 months\*

EZ Cure™ is a cross-linked collagen membrane. Guided bone regeneration promotes the proliferation of osteogenic cells in a defined area to prevent conjunctive or epithelial cells to colonize the surgical site. EZ Cure™ has the ability to maintain its functional integrity for more than 12 weeks, allowing the regenerative process to act.

#### Biocompatible and Safe

The membrane is made of porcine collagen type I and type III with standardized purification process.

EZ Cure™ collagen membrane is obtained by controlled manufacturing processes and international recognized standards for excellent biocompatibility. The collagen is extracted from veterinary certified pigs skin and is carefully purified. EZ Cure™ is sterilized by irradiation and packaged in double pouches.

#### Easy to Use Membrane

EZ Cure™ is easy to handle and conforms perfectly to the surrounding bone. Excellent handling characteristics make the procedure easiest and effective with no need to tack or suture.

### Benefits


- Biocompatible and safe
- Easy to Use
- Cells occlusive

### Indications

EZ Cure™ is intended for use in periodontal and dental surgery procedures as a material for placement in the area of periodontal defect, dental implant, bone defect or ridge reconstruction to aid in wound healing post surgery.

Considering EZ Cure™ indications and resorption time, it is recommended to combine the membrane with bone graft to new bone healing by osteoconduction (ex MBCP+™, In'Oss™)


In vitro resorption of collagen membranes study conditions, as per E-An-007 \*




More mass lost is important, more resorption in the membrane is fast. EZ Cure™ present a stable and controlled resorption.



Lateral sinus augmentation with EZ Cure™ membrane





After 3 months



After 6 months

Dr. S. Bechara, Lithuania

Resorption of the membrane EZ Cure™ in vivo \*



\* Data on files, Biomatlante

Manufacturer:  
**Biomatlante**  
 5 Rue Edouard Belin  
 ZA Les Quatre Nations  
 44360 Vigneux de Bretagne - France

[www.biomatlante.com](http://www.biomatlante.com)

ISO 13485  
 Read the instructions for use  
 Medical Device Class: III

CE 0499

Distributed by:

